

## Quality Control software for makers of immunoassay test kits.

- IMMUNQC™ is browser-based and multi-user, with a central transactional database for data storage.
- IMMUNQC™ allows for multiple user roles.
- Administrator can create an unlimited number of users.
- User account creation complies with GMP requirements.
- All user logins and system transactions resulting in data modification are logged in the Transaction
- History and are accessible to administrators through the Transaction History Log.

We'll work with you to

- Improve and implement product inserts
- Simplify your QC process
- Store, access, effectively use your historical QC data
- Train your staff to use IMMUNQC™ effectively
- IMMUNQC™ special Quality Control features include:
  - Implement package insterts; calibrators; curve, QC samples; validation rules.
  - Store reagent lot numbers, expiration dates and calibration values.
  - Maintain assay run log sheets.
  - Select lot number for assay run.
  - Allow assay run modifications to layout template.
  - Levey-Jennings chart with Westgard rules.
  - Assay run result sheet.

Create Cal Series	
ID:	
Description:	Create one or more calibrators for a particular analyte
Part Number:	AC-5701 Available!
Name:	Cal Series
Description:	
Comment:	
Raw Unit:	OD Required
Description:	Optical Density
Dose Unit:	ng/mL Required
Description:	nanograms per milliliter
Number To Create:	6
Show Items in Kit Manifest?:	<input type="checkbox"/>
Dilution Series?:	<input checked="" type="checkbox"/>
Highest Cal Amount:	100
Dilution Factor:	2.0
Include Zero?:	<input checked="" type="checkbox"/>
Conversion Factor?:	<input type="checkbox"/>